

New England Biolabs Product Specification

Product Name: dsRNA Ladder
Catalog #: N0363S
Shelf Life: 24 months
Storage Temp: -20°C
Composition: 10 mM Tris-HCl, 1 mM EDTA, (pH 8.0 @ 25°C)
Specification Version: PS-N0363S v2.0
Effective Date: 22 Mar 2021

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Marker) - The banding pattern of dsRNA Ladder on a 6% TBE gel shows discrete, clearly identifiable bands at each band of the marker, when stained with SYBR Gold at a concentration of 1 X.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 22 Mar 2021

Derek Robinson
Director, Quality Control

