
New England Biolabs Product Specification

Product Name: Low Range ssRNA Ladder
Catalog #: N0364S
Shelf Life: 24 months
Storage Temp: -20°C
Composition: 20 mM Sodium Citrate , 1 mM EDTA, (pH 6.0 @ 25°C)
Specification Version: PS-N0364S v2.0
Effective Date: 28 Aug 2019

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Marker) - The banding pattern of Low Range ssRNA Ladder on a 6% denaturing TBE-Urea gel shows discrete, clearly identifiable bands at each fragment of the marker, when stained with SYBR Gold at a concentration of 1 X.
--



Date 28 Aug 2019

Derek Robinson
Director of Quality Control

