
New England Biolabs Product Specification

Product Name: Universal miRNA Cloning Linker
Catalog #: S1315S
Shelf Life: 36 months
Storage Temp: -20°C
Composition: Supplied as a lyophilized triethylammonium salt
Specification Version: PS-S1315S v1.0
Effective Date: 06 Nov 2018

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern - The banding pattern of Universal miRNA Cloning Linker on a 15% TBE-Urea gel is evaluated against the phosphorylated starting material, resulting in a shift confirming adenylated product.



Date 06 Nov 2018

Derek Robinson
Director of Quality Control

