
New England Biolabs Product Specification

Product Name: Oligo d(T)₁₈ mRNA Primer
Catalog #: S1316S
Shelf Life: 36 months
Storage Temp: -20°C
Composition: Supplied as a lyophilized triethylammonium salt
Specification Version: PS-S1316S v1.0
Effective Date: 01 Nov 2018

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Oligonucleotide) - The banding pattern of Oligo d(T) ₁₈ mRNA Primer on a 15% TBE-Urea gel is evaluated against a control lot for expected size, sharpness, and relative intensity.
--



Date 01 Nov 2018

Derek Robinson
Director of Quality Control

