
New England Biolabs Product Specification

Product Name: Monarch[®] RNA Purification Columns
Catalog #: T2007L
Shelf Life: 24 months
Storage Temp: 25°C
Specification Version: PS-T2007L v2.0
Effective Date: 23 Mar 2021

| Assay Name/Specification (minimum release criteria) |
|---|
|---|

| |
|---|
| Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] RNA Purification Columns are tested with 10 µg of Lambda/BstEII Marker resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples. |
|---|

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 23 Mar 2021

Derek Robinson
Director, Quality Control

