
New England Biolabs Product Specification

Product Name: Monarch[®] RNA Cleanup Columns (10 µg)
Catalog #: T2037L
Shelf Life: 24 months
Storage Temp: 25°C
Specification Version: PS-T2037L v3.0
Effective Date: 28 Jun 2021

| |
|--|
| Assay Name/Specification (minimum release criteria) |
|--|

| |
|---|
| Functional Testing (RNA Recovery and Purity) - Twenty-four Monarch [®] RNA Cleanup Columns (T2037, 10 µg) are tested with 1 µg of RNA resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples. |
|---|

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 28 Jun 2021

Derek Robinson
Director, Quality Control

