

New England Biolabs Product Specification

Product Name: Monarch[®] gDNA Binding Buffer
Catalog #: T3014S/L
Shelf Life: 36 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T3014S/L v3.0
Effective Date: 21 Aug 2024

Assay Name/Specification (minimum release criteria)

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch[®] Spin Columns S2C are tested with 10 µg of HindIII digested Lambda DNA resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 21 Aug 2024

Lauren Brown
Quality Approver

