
New England Biolabs Product Specification

Product Name: Monarch[®] gDNA Wash Buffer
Catalog #: T3015S/L
Shelf Life: 24 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T3015S/L v2.0
Effective Date: 11 Jun 2021

Assay Name/Specification (minimum release criteria)

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] gDNA Purification Columns are tested with 5 µg of HindIII digested Lambda resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.
--

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 11 Jun 2021

Derek Robinson
Director, Quality Control

