
New England Biolabs Product Specification

Product Name: NEBExpress[®] *E. coli* Lysis Reagent
Catalog #: P8116S/L
Shelf Life: 24 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-P8116S/L v1.0
Effective Date: 15 Jun 2021

Assay Name/Specification (minimum release criteria)

Functional Testing (Lysis Reagent) (vGFP, Fluorescence) - A 40 µl reaction in the presence of <i>E. coli</i> cells expressing vGFP and NEBExpress [®] <i>E. coli</i> Lysis Reagent incubated for 20 minutes at 25°C results in an increase in fluorescence $\geq 7,500$ RFU in the soluble fraction compared to untreated cells.
--

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 15 Jun 2021

Derek Robinson
Director, Quality Control

